

# WIDEX

## *high definition hearing*

### Widex audio shoes

**Introducing Widex range of audio shoes in two new models for paediatric and adult fittings.**

Audio shoes can be connected to audio cables, CROS adaptors and FM receivers. The Widex audio shoe range comes with or without a lock for paediatric and adult fittings and is available for the entire Widex hearing aid range. Widex audio shoes fit firmly on to the base of the hearing aid which provides a reliable contact for constant signal transmission.

	W5	W6	W7	W8
Hearing aid type:	BV19, FL19 SV19, AK19, SD19, IN19	BV38, SV38	BV9, FL9, SV9, AK9, SD9, IN9, élan	Senso C+8/9/18/19
Audio shoe with lock	N/A			
Audio shoe without lock		N/A		

### Audio Covers

Audio covers must be used with the W6 and W8 shoes to secure the contacts on the side of the hearing aid. These do not need to be ordered separately and will be provided automatically if this audio shoe is ordered.

Audio cover (with one hole) for Senso C+ BTEs.



Audio cover (with two holes) for BV38 and SV38.

